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Referências

1. Comitê de Cuidados Respiratórios da Sociedade Torácica Chinesa. Zhonghua Jie He He Hu Xi Za Zhi. [Expert consensus on preventing nosocomial transmission during respiratory care for critically ill patients infected by 2019 novel coronavirus pneumonia]. <https://www.ncbi.nlm.nih.gov/pubmed/32077661>, 20 de fevereiro de 2020 20; 17(0):E020.
2. Organização Mundial de Saúde (OMS). Clinical management of severe acute respiratory infection when novel coronavirus (2019-nCoV) infection is suspected: Interim guidance. https://www.who.int/docs/default-source/coronaviruse/clinical-management-of-novel-cov.pdf?sfvrsn=bc7da517_2, número de referência da OMS: WHO/nCoV/Clinical/2020.3, 28 de janeiro de 2020.
3. Organização Mundial de Saúde (OMS). Infection prevention and control during health care when novel coronavirus (nCoV) infection is suspected: Interim guidance. <https://apps.who.int/iris/rest/bitstreams/1266296/retrieve>
4. Interim Clinical Guidance for Management of Patients with Confirmed Coronavirus Disease (COVID-19). Centro de Controle e Prevenção de Doenças. <https://www.cdc.gov/coronavirus/2019-ncov/hcp/clinical-guidance-management-patients.html>
5. Xiaobo Yang*, Yuan Yu*, Jiqian Xu*, Huaqing Shu* et al. Clinical course and outcomes of critically ill patients with SARS-CoV-2 pneumonia in Wuhan, China: a single-centered, retrospective, observational study. Lancet Respir [https://www.thelancet.com/journals/lanres/article/PIIS2213-2600\(20\)30079-5/fulltext](https://www.thelancet.com/journals/lanres/article/PIIS2213-2600(20)30079-5/fulltext), 24 de fevereiro de 2020. pii: S2213-2600(20)30079-5. doi: 10.1016/S2213-2600(20)30079-5.
6. Rello, J., Tejada, S., Userovici, C., Arvaniti, K., Pugin, J., Waterer, G. Coronavirus Disease 2019 (COVID-19): A critical care perspective beyond China. Anaesth Crit Care Pain Med. <https://www.sciencedirect.com/science/article/pii/S2352556820300333?via%3Dihub>, 3 de março de 2020. pii: S2352-5568(20)30033-3. doi: 10.1016/j.accpm.2020.03.001
7. Guan, W.J., Ni, Z.Y., Hu, Y., Liang, W.H. et al. Clinical Characteristics of Coronavirus Disease 2019 in China. N Engl J Med. <https://www.ncbi.nlm.nih.gov/pubmed/32109013>, 28 fev 2020. doi: 10.1056/NEJMoa2002032.

8. Josh Farkas. COVID-19. Internet Book of Critical Care (IBCC), https://emcrit.org/ibcc/covid19/#noninvasive_respiratory_support, 2 de março de 2020.
9. David S. Hui, Benny K. Chow, Thomas Lo *et al.* Exhaled air dispersion during high-flow nasal cannula therapy versus CPAP via different masks. *Eur Respir J.* <https://erj.ersjournals.com/content/53/4/1802339.long>, 2019; 53: 1802339 [https://doi.org/10.1183/13993003.02339-2018].
10. Basem M. Alraddadi, Ismael Qushmaq *et al.* Noninvasive ventilation in critically ill patients with the Middle East respiratory syndrome. *Influenza Other Respir Viruses.* <https://onlinelibrary.wiley.com/doi/epdf/10.1111/irv.12635>, 2019; 13:382–390.
11. (<https://www.usa.philips.com/healthcare/product/HCNOCTN113/respironics-niv-filter-circuit/overview>)
12. Frank Diamond. Best Approach to Disinfecting Surfaces Amid Novel Coronavirus Outbreak, *Infection Control Today*, <https://www.infectioncontroltoday.com/environmental-services/best-approach-disinfecting-surfaces-amid-novel-coronavirus-outbreak>, 10 de fevereiro de 2020.
13. Richard Lowe. Select Effective Disinfectants for Use Against the Coronavirus That Causes COVID-19. *Infection Control Today*, <https://www.infectioncontroltoday.com/environmental-services/select-effective-disinfectants-use-against-coronavirus-causes-covid-19>, 26 de fevereiro de 2020.