

Temperature Measurement




Temperature Measurement

Product #	Description	Packaging	Use with:*
1 21091A	Disposable sterile skin surface temperature probe Length: 30" (76cm) Shelf life: at least 12 months	1 box = 20 probes	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A Multi-Measurement Modules: M3000A, M3001A, M3002A SureSigns: 863064, 863065, 863066, 863068, 863087, 863088, M3015A, M3015B MP2 Monitor: M8102A, MP5 Monitor: M8105A
2 M1837A	Esophageal/Rectal Temperature Probe, 9 FR Length: 36" (91cm), sterilized, continuous monitoring Shelf life: at least 12 months		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
3 21090A	Esophageal/Rectal Temperature Probe, 12FR Length: 36" (91cm) sterilized, continuous monitoring Shelf life: at least 12 months		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement



1



2

Temperature Measurement

Disposable Esophageal/Stethoscope Temperature Probes (meets ISO 80601-2-56:2009 Standards)

Product #	Description	Packaging	Use with:*
1 21093A	Esophageal/Stethoscope Temperature Probe, 12 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months	20 probes	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
2 21094A	Esophageal/Stethoscope Temperature Probe, 18FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months		
21095A	Esophageal/Stethoscope Temperature Probe, 24 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement



Temperature Measurement

Disposable Foley Catheter Temperature Probes *(meets ISO 80601-2-56:2009 Standards)*

Product #	Description	Packaging	Use with:*
1 M2255A	Foley Catheter Temperature Probe, 14FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months	10 probes	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
2 21096A	Foley Catheter Temperature Probe, 16FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months		
3 21097A	Foley Catheter Temperature Probe, 18FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement



1



2

Temperature Measurement

Reusable Extension Adapter Cables for Disposable Temperature Probes *(meets ISO 80601-2-56:2009 Standards)*

Product #	Description	Packaging	Use with:*
1 21082A	Adapter cable: Connects disposable temperature probes to Philips (or HP or Agilent) monitors Accepts 1/8" mini phone plug. (2-pin connector) Length: 10ft (3m)	1 cable	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
21082B	Adapter cable: Connects disposable temperature probes to Philips (or HP or Agilent) monitors Accepts 1/8" mini phone plug. (2-pin connector) Length: 5ft (1.5m)		
2 989803162591	Adapter cable: Connects disposable temperature probes to Series 400 monitors Length: 10ft (3m)		Series 400 Temperature Monitors with 1/4" phone plug
2 989803162601	Adapter cable: Connects disposable temperature probes to Series 400 monitors Length 5ft (1.5m)		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement




Temperature Measurement

Reusable Autoclavable Temperature Probes for Continuous Monitoring (*meets ISO 80601-2-56:2009 Standard*)

Product #	Description	Packaging	Use with:*
1 21075A	Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 2-prong plug, reusable, continuous monitoring Length: 10ft (3m)	1 probe	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, C3: 862474, 862478, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
2 21076A	Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 2-prong plug, reusable, continuous monitoring Length: 10ft (3m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, MP2 Monitor: M8102A, MP5 Monitor: M8105A; Multi-Measurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B
3 21078A	Skin Surface Temperature Probe (autoclavable) 0.37" (9.5mm) diameter 2-prong plug, reusable tape-on, continuous monitoring Length: 10ft (3m)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Temperature Module: M1029A, SureSigns: 863064, 863065, 863066, 863068, 863086, 863087, 863088; MP2 Monitor: M8102A, MP5 Monitor: M8105A; MultiMeasurement Modules: M3000A, M3001A, M3002A, M3015A, M3015B, 860335, 860338
4 989803162621	Esophageal/Rectal Temperature Probe (autoclavable), 12FR, 1/4" phone plug, reusable, continuous monitoring Length: 10ft (3m)		Series 400 Temperature Monitors with 1/4" phone plug

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement



Temperature Measurement

Reusable Temperature Probes & Disposable Covers for SureSigns VS1/VS3/VS4 Monitors

Product #	Description	Packaging	Use with:*
1 989803162631	Esophageal/Rectal Temperature Probe (autoclavable), 10FR, 1/4" phone plug, reusable, continuous monitoring Length: 10ft (3m)	1 probe	Series 400 Temperature Monitors with 1/4-inch phone plug.
2 989803162641	Skin Surface Temperature Probe (autoclavable) 0.37" (9.5mm) diameter 1/4in phone plug, reusable tape-on, continuous monitoring Length: 10ft (3m)		
3 989803163281	Adapter cable required to connect reusable temp probes with 1/4" phone plug to legacy HP/Agilent and Philip 2-pin input connectors Length: 10ft (3m)	1 cable	Temperature probes 989803162621, 989803162631, 989803162641
4 989803163271	Adapter cable required to connect reusable temp probes with 2-pin connector to Series 400 temperature monitors with 1/4" phone plug input connector Length: 10ft (3m)		Temperature probes 21075A, 21076A, 21078A
989803179611	Genius 2 Probe Cover (CE marked), disposable, non-sterile, 2 yr min shelf life, size: 17" x 12" x 5.5"	2112 probe covers per case	VS3 (863071 #T02 & 863072 #T02)
989803179381	Genius 2 Probe Cover (not-CE marked), disposable, non-sterile, 2 yr min shelf life, size: 17" x 12" x 5.5"		VS4 (863283 #T02); 866149

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Temperature Measurement



Temperature Measurement

Reusable temperature probes & covers for SureSigns VM4/VS2/VS3/VS4 and IntelliVue MP5 monitors

Product #	Description	Packaging	Use with:*
1 989803143381	Reusable Oral Temp Probe & Well Kit predictive temp monitoring	1 probe	SureSigns VS3: 863070, 863073, 863074, VM4: 863063, VM6: 863064, 863065, VM8: 863066, IntelliVue MP5: M8105A, M8105AT, M8105AS, VS2: 863081, 863082, 863276, 863278
2 989803143391	Reusable Rectal Temp Probe & Well Kit predictive temp monitoring		SureSigns VS3: 863070, 863073, 863074, VM4: 863063, VM6: 863064, 863065, VM8: 863066, IntelliVue MP5: M8105A, M8105AT, M8105AS, VS2: 863081, 863082
M4823A	Temperature Probe Covers	1,000 covers	989803143381, 989803143391 temperature probes, 863276, 863278
989803180831	Genius 2 tympanic temperature probe, reusable, non-sterile	1 probe	VS3 (863071 #T02 & 863072 #T02) VS4 (863283 #T02); 866149

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.