

Anesthesia-Gas



Anesthesia-Gas

Product #	Description	Packaging	Use with:*
1 13901A	Reusable Nafion sample tube Length: 8 ft (2.4m)	1 tube	Module: M1026A
2 13902A	Disposable elbow airway adapter > Single patient use	20 adapters	Modules: M1026B, M1013A (G1), M1019A (G5), 866173 (G7)
13904A	Disposable bacteria filter (0.45 micron) for single patient use	20 filters	Module: M1026A
13905A	Disposable hybrid Nafion polyethylene sample tube Nafion length: 6 ft (1.8m) Polyethylene length: 9 ft (2.7m)	1 tube	
3 13907A	Calibration tube assembly	1 assembly	Modules: M1025A, M1025B, M3015A, M3015B
15210-64010	Gas: Cal 1 cylinders for tcpCO ₂	6 cylinders	M3015A, M3015B
15210-64020	Gas: Cal 2 Cylinders for tcpCO ₂		
4 M1612A	Straight airway adapter Single patient use	20 adapters	Modules: M1019A (G5), M1013A (G1), M1026A, M1026B, 866173 (G7)
M1655B	Reusable gas exhaust return line Length: 12 ft (3.6m)	1 line/5 filters	Modules: M1013A (G1), M1019A (G5), M1026B, 66173 (G7) (Note: M1026B not in the US)
M1656B	Disposable gas exhaust return filter (use with M1655B)	20 filters	Modules: M1013A (G1), M1019A (G5), M1026A, M1026B, 66173 (G7) (Note: M1026B not in the US)
5 M1657B	Water trap (use with M1658A)	25 water traps	Modules: M1013A (G1), M1019A (G5), M1026B, M1026A, 66173 (G7) with Opt. A02/A05
M1658A	Gas sample tubing for single patient use (use with M1657B)	20 tubes	
M1659A	Calibration bag (use with M1660A and M1662A)	1 bag	Modules: M1013A (G1), M1019A (G5), M1026B
M1662A	Anesthetic CAL gas Volume of gas = 6 liters (use with M1659A) Shelf life: 1 year	1 cylinder	Monitor Modules: M1013A (G1), M1019A (G5), M1026A, M1026B
6 989803191081	Philips water trap	12 water traps	866173 (G7)

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.